


pH sensor for the measurement of pH in high purity water

Swansensor pH FL

Only in combination with Swansensor Reference FL, A-87.860.100.

Technical data:

Operative and measuring range:	1 - 12 pH
Reference electrode:	Reference FL
Operating temperature:	0 50 °C
Pressure:	pressure-free
Min. conductivity:	0.055 µS/cm
Flow speed:	5 to 10 l/h
Case material:	PETP/POM
Connection:	plug PG 13,5
Weight:	40 g
Protection degree:	IP 68



Order number	Swansensor pH FL	A-87.150.200
	Swansensor Reference FL	A-87.860.100